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Published by:

USAID's Physical Rehabilitation Activity (PRA)

Handicap International, Nepal Office

233 Sallaghari Marg, Bansbari

Kathmandu, Nepal

🌐 www.hi.org 📞 +977-1-4378482 / 4374609

✉ info@nepal.hi.org 🐦 @HI_INNESKA 📺 @HI.INNESKA

Photo Sources:

USAID's PRA team

Partner Organizations - Physical Rehabilitation Centers (PRC)

Google search engine

Extracted from HI's Assistive Device Information Manual

World Health Organization Guidelines and Manuals

"This Mobility Assistive Products Service Catalogue is made possible by the generous support of the American people through the United States Agency for International Development (USAID). The contents of this catalogue are the responsibility of HI and do not necessarily reflect the views of USAID or the United States Government."

Foreword

There have been several initiatives undertaken by Nepal over the last few years to improve access to Assistive Technology (AT). These initiatives consist of developing a Nepal-specific priority list of assistive products. As a result, Nepal is the second country after Tajikistan in terms of having a priority list of assistive devices. This also includes working with the provincial and local governments in budgeting assistive products, empowering organizations of people with disabilities and other key stakeholders to build awareness and advocate for AT.

While these efforts are being made, the awareness and understanding of assistive products still present a significant challenge in the community and for those who are responsible for developing guidelines and policies for them. In addition, the diversity of products makes it difficult to fully comprehend how assistive products work. Users and relevant stakeholders need to be constantly informed about assistive technology, products, systems, and the wide range of options available to them in order to make informed decisions. In keeping with the above vision, Handicap International developed an "Assistive Device Information Manual" in 2006, which would need to be updated to reflect the current scenario.

Taking part in USAID's Physical Rehabilitation Activity (PRA), Handicap International realized that there was a need for a guiding document in the form of an assistive product catalogue for physical rehabilitation centers. This has to be used for advocacy and for raising public awareness among diverse groups of people about appropriate planning and implementation of assistive products services. Since the focus of PRA is on physical rehabilitation, the catalogue has been focused on mobility assistive devices to provide an informative and descriptive catalogue for available mobility devices. Further this publication is aimed at supporting the partners' lobbying efforts with government and non-government organizations in order to raise funds for assistive device services.

As a result of the publication of this catalogue, I am confident that this will support Physical Rehabilitation Centers (PRC) in raising awareness of mobility assistive products. I would like to congratulate all those who work in PRCs for their significant contributions towards providing quality rehabilitation services to communities. In addition, I would like to express my sincere gratitude to the authors for presenting this catalogue as a tool for the generation of uniform understanding of mobility assistive products. I would like to extend a special thank you to USAID for its generous support.

.....
Raju Palanchoke
Acting Country Director
Handicap International

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About this booklet

This booklet has been developed as an update on the existing Assistive Devices Information Manual published by Handicap International in 2006. The need has been identified for increasing awareness on mobility assistive devices to all the relevant stakeholders of physical rehabilitation service providers. The booklet includes the detail information on Assistive technology, its importance, varieties of products and information regarding mobility assistive products service providers in Nepal.

SANDIP'S STRENGTHENED DETERMINATION

© Pralhad Gairapipli / HI

“Now, I want to read and get help to improve my mobility. I need to study hard so I can get a job and become independent. I have to turn my dreams into reality! Education is my new ambition! I want to study hard and qualify for jobs. I plan to open a mobile repair shop or start working after I complete my education.”

PROSTHESIS SHOULD BE ADAPTED AS CHILD GROWS

© CBRB / HI

"Growth is an important aspect of a child's life and as their bodies change, prostheses have to be adapted or changed in the similar manner to accommodate them. Just as they outgrow shirts, pants, and shoes, they will outgrow their prostheses."

List of Acronym

AFO	Ankle Foot Orthosis
AT	Assistive Technology
EDCD	Epidemiology and Disease Control Division
EVA	Ethylene Vinyl Acetate
FO	Foot Orthosis
GATE	Global Cooperation on Assistive Technology
ICRC	International Committee of Red Cross
KAFO	Knee Ankle Foot Orthosis
KO	Knee Orthosis
LCDMS	Leprosy Control and Disability Management Section
MOHP	Ministry of Health and Population
MOWCSC	Ministry of Women, Children and Senior Citizen
PAPL	Priority Assistive Product List
PRA	Physical Rehabilitation Activity
PRC	Physical Rehabilitation Center
UHC	Universal Health Coverage
USAID	United States Agency for International Development

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Assistive technology can be defined as “any piece of equipment, or product, whether it is acquired commercially, modified, or customized, that is used to increase, maintain, or improve the functional capabilities of individuals with disabilities”. Assistive technology ranges from mobility devices such as wheelchairs, crutches, prosthesis, orthosis, etc. to hearing aids, white cane or special devices to support activity of daily living like dressing, feeding, washing and more).

Improving access to high-quality and accessible assistive products and equipment, for which five areas are interdependent on each other.

- **People:** user-centric; Work closely with users and user groups
- **Policy:** Assistive Technology Policy Development
- **Products:** Development of National Primary assistive product List
- **Provision:** Flow system and model suitable service guideline
- **Personnel:** Training and Skilled workforce Development

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Mobility assistive devices are subset of AT that are designed to facilitate or enhance a user's personal mobility – this relates to their ability to change and maintain body position and walk and move from one place to another. Mobility devices are one of the most common types of assistive technologies, products or devices.

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IMPORTANCE OF MOBILITY ASSISTIVE DEVICES

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Mobility assistive devices are important

- To have a significant impact on the level of independence and participation which people with disabilities are able to achieve
 - To reduce the need for formal support services as well as reduce the time and physical burden for caregivers.
 - To enable individuals to lead healthy, productive, independent and dignified lives and to participate in various aspect of their life such as education, livelihood and social aspects
-

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SERVICE USERS OF MOBILITY ASSISTIVE DEVICES

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People that can benefit from mobility devices are those who experience mobility difficulties as a result of a broad range of health conditions and impairments, including amputation, arthritis, cerebral palsy, poliomyelitis, muscular dystrophy, spinal-cord injury, spina bifida, stroke and visual impairment.

These device are are also relevant for elderly people who experience mobility difficulties.

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

GATE Global Cooperation on Assistive Technology

Only 1 in 10 people in need has access to assistive products.

0.5% population in developing countries are in need of prostheses/orthoses and 1% require wheelchair.

Initiation of Global Cooperation on Assistive Technology (GATE) as a global initiative to assist Member States to improve access to assistive technology as a part of Universal Health Coverage. Some international documents are shown here.

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Appropriate technology is a system that provides proper fit and alignment based on sound biomechanical principles, which suit the needs of the individual and can be sustained by the country at the most economical and affordable price.

An Appropriate technology for Mobility Assistive Devices must

- Meet the user's needs and environmental conditions;
- Provide proper fit and support;
- Be safe and durable;
- Be available in the country: and
- Be obtained and maintained and services sustained in the country at an affordable cost.

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
Pricing is an essential step and should be done in a scientific way. Therefore, when pricing assistive devices, the following points are taken into consideration:

- Proper review of all direct and indirect costs
- Use of transparent and scientific method
- The price varies as there are many types of products. The more complex the product, the higher the cost
- As production, delivery processes and workforce play an important role in assistive product services, in depth analysis is done.
- As raw materials and components are imported from abroad, it is affected by national and international policies and regulations.

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Inappropriate prescription or provision of assistive products lead to:

- Discomfort, pain
 - Risk of injuries, wounds and pressure sores
 - Secondary complications such as contractures, deformities, etc.
 - Abandonment of the devices due to improper fit or adjusted environmentally
-

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ASSISTIVE DEVICE SERVICE PROVISION

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
Service delivery steps

Assistive devices service provision must follow the above shown four steps.

Assessment: Should be assessed by a professional with appropriate knowledge and skills or, when possible and required, by a multidisciplinary team. Rehab Doctor, Prosthetist and Orthotist, Physiotherapist, Occupational therapist, Social Worker, etc.

Fabrication and Fitting: Ensure that the full potential of products is used and hazards for users and personnel are minimized

User Training: Training is critical to the outcome of treatment and helps to maximize the benefits of products

Product delivery and Follow-up: Before delivery, all features must be checked according to a protocol to ensure the best possible quality, appropriateness and function. Follow-up is required at certain intervals to review outcomes and ensure that there are no problems with the fit, comfort, function or durability of the device.

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Some assistive products need to be assembled (put together), fabricated (made) or adjusted for the user whereas some are available prefabricated or need to be unboxed and assembled. Trained personnel prepare the products for an initial fitting.

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RECENT DEVELOPMENT ON ASSISTIVE TECHNOLOGY IN NEPAL

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Rehabilitation and Assistive Technology are included in policy documents of Nepal such as Nepal Health Sector Strategy, 10 years strategy and action plan for disability management, Priority Assistive Product List (PAPL) of Nepal, etc.

Leprosy Control and Disability Management Section (LCDMS) is the designated focal point for Disability and Rehabilitation including Assistive Technology for Ministry of Health and Population.

Publish of several guidelines, standards and protocols is key step towards integration.


Some funds allocated by government of Nepal through LCDMS, Ministry of Women, Children and Senior Citizen (MOWCSC), provincial and local authorities for assistive product services.

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COMMON TYPES OF MOBILITY ASSISTIVE DEVICES

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Priority Assistive Product List (PAPL) of Nepal has already categorized mobility assistive products into

- Walking Aids: Crutches, Walking Stick, Walkers
- Wheelchairs: Manual/Electric/Sports Wheelchairs, Tricycle
- Prosthesis : Lower limb prosthesis, Upper limb prosthesis
- Orthosis: Lower Limb Orthosis, Upper Limb Orthosis, Spinal Orthosis
- Developmental Aids: Special chair, Standing Frame, Prone board

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
Photo credit @ physiopedia

Prostheses are the external device which replaces the whole limb or missing segment, i.e. given for functional, replacement of missing segment or cosmetic reasons.

Beneficiaries

People who have completely or partly lost their arms/legs either by birth or after birth irrespective of the causes

Indications



Abnormal upper or lower limb (by birth); Amputations due to accident, diabetes, leprosy, frostbite etc.

Considerations

Special care to be taken for people having issues with blood circulation (e.g. diabetes) or with sensation (e.g. leprosy)

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This prosthesis is prescribed for people with partial foot amputation or missing part of the foot. In partial foot amputation the part of the foot is removed due to any particular condition may it be trauma, disease, etc.

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This is developed based on individual requirements.
Can be made up of laminated or thermoplastics.

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It is prefabricated and made up of
Silicone material.

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
Toe fillers are used as an insert in the footwear and used
for the replacement of toes. It also improves the life span of
footwear and provide support.

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This prosthesis is prescribed for the people with amputation/missing through the ankle joint.

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
This is developed based on individual requirements. Due to bulbous nature of the stump one design of syme's prosthesis is made with medial window or posterior opening for easy put on and off. The other design is with push fit. Can be made up of laminated or thermoplastics. Specific low profile foot are used for length management.

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This prosthesis is prescribed for the people with below knee amputation or missing part of the leg below the knee joint.

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It contains Socket, shank and foot. It simulate the natural structure of human limb as weight is transferred through central structural segment. The components are purchased and are of various types while socket is locally developed as per individual measurements.

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
The socket and shank is made up of same material and do not have central structural segment and weight is transferred through surface. The prosthesis is locally developed as per individual measurements. Various types of foot can be used.

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This prosthesis is prescribed for the people with through knee amputation or missing part of the leg through the knee joint.

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It contains Socket, knee joint, shank and foot. The components are purchased and are of various types while socket is locally developed as per individual measurements.

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The socket and leg section is connected through external hinge and bars and do not have central structural segment and weight is transferred through surface. The prosthesis is locally developed as per individual measurements. Various types of foot can be used.

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This prosthesis is prescribed for the people whose all segment below the thigh is removed or missed. For this amputation, the prosthesis must replace the knee, leg and foot segment.

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This prosthesis is prescribed for the people whose all segment below the hip joint is removed or missed. For this amputation, the prosthesis must replace hip joint, thigh, knee joint, leg and foot segment.

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It contains Socket, hip joint, thigh, knee joint, shank and foot. The components are purchased and are of various types while socket is locally developed as per individual measurements.

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Used for a person who has an amputation through or close to the ankle. The heel part of the foot does not come up as high as other feet so there is more room to fit prosthetic components without making the leg too long

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A basic and affordable foot. There is no movement at the ankle. This reduces moving parts (making maintenance easier), however can make it difficult for the user to walk over uneven ground or up and down slopes.

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Provides some ankle movement, providing more shock absorption. The user can walk with a smoother walking pattern, and more easily over uneven ground or up and down slopes.

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
Gives the user more ankle movement than the single axis foot.

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Provides more function for the user. Provides good shock absorption, stability over uneven ground and acts like a spring to push the foot off the ground for a more natural walking pattern

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A basic knee joint. Some basic single axis knees will be able to be locked in a straight position for extra stability, however this can make normal walking more difficult as the knee does not bend.

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
More adjustability than a basic single axis knee. The amount of stability and how fast the knee bends can be adjusted to suit the user. It can also help to straighten the knee during walking. It can be mechanical, pneumatic or hydraulic controlled. Several designs and functions are available.

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Works more like a real knee. It looks better when sitting, and is more stable when walking than other knees. This knee is good for people who have an amputation through the knee, very short stumps or have weak hip muscles. Several designs and functions are available.

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
Single axis knee with manual locking system. Used with ICRC modular system.

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
This prosthesis is prescribed for the people with partial hand amputation or missing part of the hand or finger. In partial hand amputation the part of the hand is removed due to any particular condition may it be trauma, disease, etc.

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This is a prefabricated product and procured. This hand is made up of Silicone material. Available in sizes, color and gender based.

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
This is a prefabricated product and procured. This finger is made up of Silicone material.

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This prosthesis is prescribed for the people with amputation/missing parts below the elbow .

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It consists of socket and functional hand segment operated by body power through attached cable control. The socket is fabricated locally through lamination or thermoplastic method whereas the below elbow kit is purchased. Hook may also be used.

सौन्दर्य (कस्मेटिक) कुहिनो देखि तलको कृत्रिम हात / AESTHETIC (PASSIVE) BELOW ELBOW PROSTHESIS

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
Here the passive hand is attached with the customized socket. It is used for aesthetic purpose.

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This prosthesis is prescribed for the people whose all segment below the arm is removed or missed. For this amputation, the prosthesis must replace the elbow, forearm and hand segment.

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It consists of socket and functional elbow and hand segment operated by body power through attached cable control. The socket is fabricated locally through lamination or thermoplastic method where as the above elbow kit is purchased. Hook may also be used.

सौन्दर्य (कस्मेटिक) कृत्रिमो देखि माथिको कृत्रिम हात/ **AESTHETIC (PASSIVE) ABOVE ELBOW PROSTHESIS**

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Orthoprosthesis is an external device which is neither a prosthesis or orthosis, but made by the combination of both prosthetic and orthotic principles

Beneficiaries

People who have significant limb length discrepancy (lower limb) and surgical intervention for amputation is not performed

Indications

Leg Length discrepancy which cannot be accommodated with shoe compensation; lower leg deformity with a “short leg” as a result

Causes of very short leg

Congenital (such as absence of bones or portion of leg or whole leg) ; Acquired (Trauma, untreated malunited fracture)

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It contains Socket/Shell, shank and foot. The discrepancy is so high that the extension has to be done by attaching prosthetic components. The components are purchased and are of various types while socket is locally developed as per individual measurements.

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The weight is transferred through surface.

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Orthosis is an external device that supports the weakened or paralysed part of the limb or whole limb

Beneficiaries

Persons who have limitations in mobility and activity of daily livings due to muscular or skeletal problems, and paralysis of limbs.


Purpose

Provide support; Prevent/Reduce deformities; Enhance Function

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
Commonly known as FO is prescribed in case of weakness and deformities of foot. Includes customized supports, insoles, orthopedic/molded shoe and components.

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This is made with leather or other shoe materials. It may be combined with several modifications such as medial arch support, raise, metatarsal bars, etc. as per condition.

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This type of foot orthosis is made up of hard material with soft lining. It is directly inserted inside the shoes. It can be either custom made or readymade.

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

This is made up of soft EVA or silicone with gel and is given for comfort.

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Commonly known as AFO is prescribed in case of weakness and deformities of lower leg, ankle and foot. AFO are of several designs as per the conditions and requirement to control the foot and ankle structure.

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
This is made with plastics and have numerous designs depending upon the need of control and mobility. Fixed AFO, Articulated AFO, Dynamic AFO, anterior, etc.

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It is plastic made, flexible in nature and prefabricated. It's popular as Foot Drop Splint .

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
Provided for children with club foot. Designs vary and can be adjustable.

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Commonly known as KO is prescribed in case of weakness and deformities related to the knee and also often utilized temporarily after different types of knee surgeries.

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This is made with plastics and have numerous designs some are with knee joint and some are fixed depending upon the need of control and mobility. It is custom made as well as readymade. It is provided for pain reduction, ligament injuries, prophylactic, knee deformities, etc.

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It is made with fabric including metal/plastic stripes for reinforcements. It is custom made as well as readymade. It is used as knee immobilizer to prevent contractures , heal injuries faster, post surgical, etc.

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
Commonly known as KAFO is prescribed in case of weakness and deformities in the upper leg and the knee, and possibly also ankle and foot.

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This is made with plastics and have numerous designs some are with knee joint and some are fixed depending upon the need of control and mobility. It is custom made. It is provided for weight bearing and mobility.

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It is made up of leather and metal and include orthopedic boot. It is custom made.

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Upper limb orthosis serve one or more basic purpose:

- To assist weakened muscle or to substitute the total loss of muscle power
- To protect the part from pain or potential deformity
- To correct existing deformity

Beneficiaries

Persons having problems with the function of their hand by birth, due to injury or any other reason, hampering their participation in daily activities.

Types


Finger Orthosis
Wrist Hand Orthosis
Elbow Orthosis

Indications

Nerve Injury
Hemiplegia
Burn Injury
Cerebral Palsy

Considerations

Special care to be taken for diabetic case

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
This is made with plastics and keeps the hand in desired static position. It is custom made or readymade.

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It is made up of plastic and wrist can be progressively positioned through the attached turn buckle.

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It is made up of plastic and dynamically positions the wrist as well as fingers into the functional position. It is also used for therapeutic purpose.

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This is made with plastics and keeps the elbow in desired static position. It is custom made or readymade. It can be jointed or fixed both or can be used with turnbuckle.

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It is made with fabric including metal/plastic stripes for reinforcements. It is custom made as well as readymade. It is used as elbow immobilizer to prevent contractures, heal injuries faster, post surgical, etc.

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BENEFICIARIES


Persons with low back pain and those who have visible deformities of back bone (Spine) either correctable or fixed.

INDICATION

Scoliosis (side wise curvature of spine);
Kyphosis (frontward curvature); Cerebral palsy; Low back pain;

CONSIDERATIONS

Special care are to be taken for muscular dystrophy and non-operative treatment of Arthritis of joint of spine. Further attention should be given in appropriate prescription and production of the orthosis to prevent complications.

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This is made with plastics and keeps the whole trunk in appropriate position preventing further progression of problems. It is custom made orthosis. This is just a design among many available which is based on the specific conditions.

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It is available in both flexible as well as rigid design. Mostly readymade products are available, however custom made can also be prepared.

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It is made up of metallic frame with flexible fabric corsets. Commonly available are Taylor, Knight Taylor , etc. Lumbar belt is flexible .

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BENEFICIARIES

Having poor balance while walking/standing; cannot or may not put all their weight on the lower limb; use calipers or prostheses on both legs

INDICATIONS

Fracture (for temporary use); Joint instability; Poor Balance; Post surgical ambulation (temporary use); person with bilateral prostheses/orthoses; Degenerative conditions accompanied by pain; Neurological problems; Deformities

CONSIDERATION

People having balancing issues; Risk of deterioration while standing or walking; Very weak or painful upper limbs; Poor coordination skills; Fresh lower limb fractures

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BENEFICIARIES

Person who cannot stand or walk for whatever reason, including: some people with prostheses/orthoses; children who are floppy and slow in developing the ability to sit
Who cannot control the movement of the body
Weak and get tired easily

INDICATIONS

People with multiple disabilities that restrict independence in mobility
People with paralysis
People with spinal cord injury
Elderly people

CONSIDERATION

Special care is needed for
People with epilepsy
Special attention to pressure sores and loss of sensation

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The user has to propel manually or with the support from the helper. There are several designs and models in manual wheelchairs ranging from indoor to outdoor mobility. Available in folding and non-folding designs with multiple sizes.

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It is rough terrain wheelchair available in several sizes. This is an active wheelchair.

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This is specially designed to provide postural support to users with higher needs in appropriate seating. Suitable for children who requires support to maintain posture. Available in sizes.

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The sport wheelchairs suit the specific needs for the sport. are usually non-folding, with a pronounced negative camber for the wheels, and often are made of composite, lightweight materials.

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It is battery operated wheelchair propelled by electric motor.

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
A foam cushion designed to provide some cushioning and comfort for a wheelchair user who is not at risk of developing a pressure sore and does not need postural support from their cushion.

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A cushion designed and proven to provide pressure relief for the prevention of pressure sores for people who are at risk of developing a pressure sore. There are different types of pressure relief cushions available including contoured foam, air and fluid filled cushions. Pressure relief cushion can be locally made using chip foam cushion as well.

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The tricycle is for outdoor mobility especially for travelling long route. It is manually operated with the help of hand pedal and steered through the handle.

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
The ramps represent one of the main equipments for universal accessibility. Ramps are the slope to enable the persons with mobility limitation to access the infrastructure.

Ramps and steps need to be of appropriate and precise size to fully assist in the physical access of people with functional limitations.

The ramps should be built according to minimum requirements and standards for safety and usability (slope, minimum width, maximum length, landing and resting zone, safety features, design...) - these technical requirements should be sought before construction.

It helps the wheelchair and other assistive device users to move up and down on their own.

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

Many children have problems with sitting/standing balance and coordination of their different daily activities. A developmental aids facilitates to develop their Gross motor functions; hand functions; communication and interaction

BENEFICIARIES

Have poor balance while sitting/standing
Have coordination problems for daily activities


INDICATIONS

Poor sitting balance
Problems to stand without support
Poor control of head
Coordination problems

CONSIDERATION

These aids may be of adjustable or fixed heights and can be fabricated according to individual needs

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Supporting proper sitting position also help children who have difficulties with their upper limb to better control them (a Cerebral Palsy child who is in a proper and stable sitting position will be able to better use his/her hands). There are also advantage in terms of interaction with other - being in a proper seated position allows to better see and look around and interact with others.

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Most children with balance or control problems for standing can play or stand with the help of a standing frame. This helps to keep the weight of the body in the legs, which helps to increase blood circulation and keep the bones strong. It is commonly made up of wood which consists of a frame to stand, lap table to play or do hand activities and back support. Small castor wheels for transfer can also be attached.

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
The prone board helps children with cerebral palsy to reduce muscle spasms, prevent contractures, and help to gain head control, as well as develop shoulder stability. It is either static or adjustable

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S.No.	Service Center Name	Address & Contact No.	Province
1	Community Based Rehabilitation, Biratnagar (CBRB)	Ishwar Marg, Biratnagar 021-531594	1
2	BP Koirala Institute of Health Sciences	Dharan, Sunsari 025-525555	1
3	Prerana Rehabilitation Center	Malangawa, Sarlahi 046-520442	2
4	Prerana Satellite Unit	Navalpur, Sarlahi 9803731301	2
5	Lalgadh Hospital	Lalgadh, Dhanusha 041-620182	2
6	Spinal Injury Rehabilitation Center,	Sanga, Kabhre 011-660847, 660848	3
7	Hospital and Rehabilitation center for disabled children	Janagal, Kabhre 011-661666, 661888	3
8	Physical Rehabilitation Center, National Disabled Fund	Bhrikuti Mandap, Kathmandu 01-4224968, 4239586	3
9	Nepali Army Rehabilitation Center	Sano Bharyang, Kathmandu 01-4880684	3
10	Limb Care Pvt. Ltd.	Bus Station, Budhanilkantha 9841301073	3
11	Orthocare Rehabilitation Center	Ekikaran Road, Kathmandu 9841267004	3
12	Mobility Nepal Pvt. Ltd.	Khusibun Bishnumati Right Side, Kathmandu 01-4389410	
13	Anandban Hospital, the Leprosy Mission Nepal,	Anandaban Hospital, Lalitpur 01-6218398 Satdobato Clinic, Lalitpur 01-5151371	3
14	Orthopedica	Kupondole, Lalitpur 01-5180349	3
15	Green Pastures Hospital	Hariyo Kharka, Pokhara 061-431162, 430342	4
16	Nepalgunj Medical College Rehabilitation Center	Kohalpur Banke 081-521572	5
17	Nepal National Social Welfare Association	Mahendranagar, Kanchanpur 099-522182	7